

DEMONSTRATION VIDEOS ON:
WWW.DERMASCULPT.NET




BLUNT-TIPPED FLEXIBLE ATRAUMATIC MICROCANNULA



DERMA **SCULPT**
BLUNT-TIPPED MICROCANNULA






CUTS THROUGH
BLOOD VESSELS
& NERVES CAUSING
BRUISING & PAIN

SHARP NEEDLE

VS

DERMASCULPT

SLIDES AROUND
VESSELS & NERVES
AVOIDING
BRUISING & PAIN



TESTIMONIALS:

ADAM MAUER MD - Garden City, NY

Thank you ! Your product changed my whole practice. Such an amazing device, what took the world so long and why was I slow in discovering it? Well, better late than never.

KIAN KARIMI MD - Santa Monica, CA

I used the DermaSculpt Microcannula to perform lip augmentation, and I was very impressed with the fact that there was nearly no bruising after the procedure [...]

JOEL L. COHEN, MD - Denver, CO

I have found the DermaSculpt cannulas to be incredibly helpful in minimizing bruising and helping to achieve outstanding results in soft tissue augmentation of the cheeks, infra-orbital hollows, dorsal hands as well as the décolletage. Once you start using these cannulas, you totally get hooked. [...]

» more testimonials on dermasculpt.net...


DEMONSTRATION VIDEOS ON WWW.DERMASCULPT.NET



COMPARED TO USING A SHARP NEEDLE

- ✓ Virtually no bruising
- ✓ Single point of entry for wider coverage
- ✓ Considerably less pain
- ✓ Production of new collagen
- ✓ Minimal downtime

**FULL
COVERAGE
W/ FEW
ENTRY
POINTS**



For medium to thick HA and non HA fillers, median to deep injections



25G x 1 1/2" (38mm)

can be used for most facial indications



25G x 2" (50mm):

extra length to cover larger treatment area

For viscous & highly concentrated fillers, injections deep into the skin to create volume



22G x 2" (50mm)

For facial fat graft (not recommended for fillers)



18G x 2 3/4" (70mm)

AVAILABLE SIZES:

For low concentration and viscosity HA fillers, superficial injections



30G x 1" (25mm):

tear trough, fine lines - precise work

For fillers of medium concentration and viscosity HA fillers



27G x 1" (25mm):

precise work



27G x 1 1/2" (38mm):

most versatile size, used for most facial indications



27G x 2" (50mm):

extra length to cover larger treatment area

MAIN ADVANTAGES OF DERMASCULPT®

- * Much safer than using a classical sharp needle because DermaSculpt® does not cut through vessels and tissues but rather navigates smoothly through the dermis and Sub-Q layer (depending on thickness and flexibility of microcannula used). The blunt end and flexibility of the microcannula allow for virtually painless and atraumatic procedures.
- * When using the fanning technique and thanks to the length of the microcannula, very few entry points are needed in order to cover the entire face.
- * There are also strong reasons to believe that the mechanical action of the DermaSculpt® microcannula on the fibroblasts generates new collagen into the skin.

1 Puncture a hole into the skin using a larger-gauge sharp needle

The gauge of the sharp needle needs to be slightly superior to the gauge of the microcannula in use.

For instance, 26G Sharp needle for 27G microcannula; 24G Sharp needle for 25G microcannula; etc

* These sharp needles are supplied together with the DermaSculpt Microcannulae.

2 Insert the tip of the microcannula into the skin perpendicularly.

Insertion of the tip of the microcannula should be from a 60 to 90 degree angle until the correct depth has been found.

Pinching the skin will create a small drop of blood out of the puncture hole which should help you locate the hole.

* Place the index finger of the injecting hand close to the tip of the microcannula to help guide the cannula into the skin. Stretch the skin opposite direction with the non-injecting hand to help enlarge the puncture hole.

THE 4 EASY STEPS OF THE SKIN SCULPTING TECHNIQUE™ (SST)

3 Slide the Dermasculpt microcannula into the skin.

As you are entering the skin, rotate the cannula tangentially to the skin surface using the flexibility of the microcannula.

If any resistance is encountered, go back and forth without forcing.

* Depending on the gauge you've selected for a specific indication or for a specific filler, the microcannula will start navigating either in the deep dermis or the sub-Q layer.

4 Inject filler in a retrograde fashion, fanning the entire area to be treated.

Start injecting while pulling the length of the microcannula in and out of the skin. To change direction, pull the cannula out of the skin almost entirely, back to step two (leaving the tip at the desired depth).

DermaSculpt is a brand of:



Aesthetic Experts Laboratory
London W1B 3HH

✉ info@aelaboratory.com
🌐 www.aelaboratory.com

Distributed by :

