

# FarrowWrap® Footpiece Application Techniques

Applicable for the following fabrics


- JOBST® FarrowWrap® Lite
- JOBST® FarrowWrap® Strong
- JOBST® FarrowWrap® Classic



1 Apply a liner. Either the open toe Soft Liner (white) or the closed toe Silver Liner (beige). With the Soft Liner, ensure you extend the fabric past the toes distally and past the knee proximally to allow sufficient fabric to fold back over the wraps once they have been applied. Unwrap the footpiece and lay it on the floor with the hole for the heel closest to you.



2 Place the foot on top of the footpiece, resting the heel over the hole in the footpiece.



3 To secure the foot, hold the footpiece flap closest to the toes on each side and pull both sides over the top of the foot to the opposite side. Apply at or near end stretch with the VELCRO®. Now apply the remaining flap (the one closest to the ankle) in the same manner.



4 After securing the forefoot, hold one side of the ankle band and pull it over the top of the mid-foot to the opposite side and secure it in place at or near end stretch.



5 Repeat the same step with the other half of the ankle band on the opposite side. When removing the footpiece, fold the VELCRO® over on itself and back on the same band.



6 **Care instructions (Strong and Lite fabric):** Machine wash in cold water with a mild detergent, tumble dry with no heat, do not iron, do not dry clean.

**Care instructions (Classic fabric):** Hand wash in cold water with mild detergent, lay flat to dry, do not iron, do not dry clean.

# FarrowWrap® Legpiece Application Techniques\*

Applicable for the following fabrics

- JOBST® FarrowWrap® Lite
- JOBST® FarrowWrap® Strong
- JOBST® FarrowWrap® Classic



1 Apply a JOBST® Farrow liner and footpiece if applicable.

Unwrap the legpiece and lay the wrap on the floor behind your leg with the label side furthest back and facing up. The shortest bands will be closer to the foot.



2 Place the wrap with its spine positioned behind the leg starting from the ankle and finishing just below the knee.

Apply the top wrap as an anchor to hold the wrap in place. Then commence wrapping from the ankle up towards the knee.



3 Grab each side of the band closest to your ankle and pull forward wrapping the band around the lowest part of the leg at or near end stretch and affix with the VELCRO®. Undo the top band before applying the second last band to ensure the final band sits over the top. If wearing a footpiece, the bottom band should overlap the footpiece.



4 Repeat with each band moving up the leg, smoothing out any wrinkles as you go.

Try to ensure a 50% overlap between adjacent straps.

**Care instructions (Strong and Lite fabric):** Machine wash in cold water with a mild detergent, tumble dry with no heat, do not iron, do not dry clean.

**Care instructions (Classic fabric):** Hand wash in cold water with mild detergent, lay flat to dry, do not iron, do not dry clean.

\* Please refer to Instructions for use, in product pack.