






Blade Range / Part No.

Blade shape	Blade dia. mm				
Circular	64	19610*	19640	19660	20240
Segment	64	19670	19680	19730	20280 *
Segment	80	19700	19720	19710	20290 *

	PTFE	TIN	ZrN	<i>NEW</i> IonN
POP - Gypsona	Optimal			
Polyester Deltalite/Benicast/Cellacast	Suitable			
Fibreglass Deltalite/Scotchcast	Not recommended			

NEW

Ion Nitride Blade range
for Fibreglass users



* Longer 80 dia. mm
segment blade, useful for
standard or **SPICA** cases

* 19610 & 20280 blade supplied with each cast saw

Accessories



Extractor Filter
(single item)
Order code 670113




CC7 hood
Order code 670713



Extractor HEPA dust bag
(pack of 10)
Order code 670123



2m Hose kit
Order code 19600



CC7 vacuum attachment
(with extractor system)
Order code 670733



Cast saw wall
mount for CC8
Order code 19850

CLEANCAST™

ADVANCED CAST SAW AND DUST EXTRACTION SYSTEM

For further information on De Soutter Medical products, please contact the address in your area.

United Kingdom (International Head Office)

De Soutter Medical Limited
Halton Brook Business Park, Weston Road
Aston Clinton, Aylesbury, Bucks HP22 5WF
Tel: +44 (0) 1296 634 000
Fax: +44 (0) 1296 634 033
Email: info@de-soutter.com

Australia

De Soutter Medical Australia Pty Ltd
2/12-14 Apollo Drive, Hallam, Victoria 3803
Tel: +61 (0) 3 9702 4441
Fax: +61 (0) 3 9702 4484
Email: australia@de-soutter.com

Austria

De Soutter Medical Austria
Zweigniederlassung Österreich
Dietrichsteingasse 10, A-3400 Klosterneuburg
Tel: +43 (0) 676 96 71 770
Fax: +43 (0) 2243 21 656
Email: austria@de-soutter.com

Belgium

De Soutter Medical Belgium
Bessemerstraat 14, 3620 Lanaken
Tel: +32 (0) 89/47 15 37
Fax: +32 (0) 89/70 12 19
Email: belgium@de-soutter.com

France

De Soutter Medical France
1 252 Avenue Parc des Expositions
33260 La Teste de Buch
Tel: +33 (0) 5 56 54 89 36
Fax: +33 (0) 9 70 61 37 60
Email: france@de-soutter.com

deSoutter
MEDICAL

Germany

De Soutter Medical Germany
Bahnhofstraße 4,
D-66625 Nohfelden
Tel: +49 (0) 68 52-99 12 46
Fax: +49 (0) 68 52-99 12 47
Email: deutschland@de-soutter.com

Italy

De Soutter Medical italia
Localita' Fornace SNC
27022 Casorate Primo - PV, Italia
Tel: +39 (0) 2 9009 4098
Fax: +39 (0) 2 9009 2673
Email: italy@de-soutter.com

Netherlands

De Soutter Medical Nederland
Gelderlandhaven 2X
3433 PG Nieuwegein
Tel: +31 (0) 85 0491480
Fax: +31 (0) 85 0491489
Email: nederland@de-soutter.com

United States of America

De Soutter Medical USA, Inc.
224 Rolling Hill Road, Suite 12A
 Mooresville, NC 28117
Tel: +1 (704) 655 9040
Fax: +1 (704) 987 2035
Email: usa@de-soutter.com