


# Cutimed® Siltec®

Closing wounds. Together.

**Super  
Absorbent  
silicone foam  
dressing**

Cutimed® Siltec® -  
The foam dressing  
that adapts to  
various exudate  
levels



**THERAPIES. HAND IN HAND.**


**BSN** medical

an SCA company

# DESIGNED TO WORK FOR YOU

## THE NEW CUTIMED® SILTEC® COMPONENTS

The right combination of specially designed materials makes Cutimed® Siltec® deliver the performance you expect from a foam dressing.




## CUTIMED® SILTEC® RANGE

### CUTIMED® SILTEC®

When it comes to choosing a dressing that can be changed without damaging the wound or skin around it, Cutimed® Siltec® is driven by technology that meets the challenge. It is a soft and highly conformable foam dressing that adapts to all wound condition.

- Low to moderately exuding wounds
- Maintains critical moisture balance
- Skin protection with low adherence while reducing risk for maceration




Cutimed® Siltec®

### CUTIMED® SILTEC® L

Cutimed® Siltec® L is a thinner version of Cutimed® Siltec® and is a soft and highly conformable foam dressing for low to moderate exuding wounds. Cutimed® Siltec® L can be changed without damage to the wound or surrounding skin.

- Soft and conformable wound dressing
- Protective dressing for compromised and fragile skin



Cutimed® Siltec® L

## CUTIMED® SILTEC® B

Cutimed® Siltec® B has been designed to facilitate gentle adherence with effective exudate management for easy dressing application and removal.



Cutimed® Siltec® B

## CUTIMED® SILTEC® SACRUM

In order to ensure reliable exudate management when treating sacral ulcers, BSN medical has developed Cutimed® Siltec® Sacrum. This bordered foam dressing is particularly shaped to fit the awkward sacral area, offering an extra-thick wound pad for a high absorption capacity, plus providing additional cushioning.

- 5mm wound pad (with Super Absorbers) for moderate exudate levels and cushioning effect
- Wide, all-round adhering border
- Silicone wound contact layer for pain free and atraumatic removal
- Particularly shaped to fit the sacral region
- Shape allows for secure adherence with reduced friction and shear when immobile patient is repositioned
- Available in 2 sizes



Cutimed® Siltec® Sacrum

## CUTIMED® SILTEC® HEEL

The dressing is designed to fit the anatomical shape of the heel and is folded around the heel at application.

- Super Absorbent Soft Silicone Foam Dressing
- Reduces pressure peaks, provides cushioning



Cutimed® Siltec® Heel 2D

## CUTIMED® SILTEC® HEEL 3D

Cutimed® Siltec® Heel 3D is especially designed to provide an optimal fit to every heel shape during ulcer management. Its anatomical shape provides cushioning and pressure relief from the heel to the Achilles tendon.

The extra thick foam of the dressing allows for reliable absorption of exudate.







- Super Absorbent Soft Silicone Foam Dressing
- Ready-to-use, easy application



Cutimed® Siltec® Heel 3D

# PRODUCT OVERVIEW

## Order information

Cutimed® Siltec®	Pad size	Wound pad	Dressing size	Ref.-no.	Dressings/box
<b>Cutimed® Siltec® L (Light)</b>  Thin silicone foam dressing with Super Absorbers			5 cm x 6 cm	7328300	10
			10 cm x 10 cm	7328301	10
			15 cm x 15 cm	7328302	10
<b>Cutimed® Siltec®</b>  Silicone foam dressing with Super Absorbers			5 cm x 6 cm	7328500	10
			10 cm x 10 cm	7328501	10
			10 cm x 20 cm	7328502	10
			15 cm x 15 cm	7328503	10
<b>Cutimed® Siltec® B (Border)</b>  Self-adherent silicone bordered foam dressing with Super Absorbers	(4 x 4 cm)		7.5 cm x 7.5 cm	7328400	10
	(6 x 6 cm)		10 cm x 10 cm	7328405	10
	(8.2 x 8.2 cm)		12.5 cm x 12.5 cm	7328401	10
	(6 x 16.5 cm)		10 cm x 22.5 cm	7328406	10
<b>The anatomical specialists</b>		<b>Wound pad</b>	<b>Dressing size</b>	<b>Ref.-no.</b>	<b>Dressings/box</b>
<b>Cutimed® Siltec® B Oval</b>  Self-adherent silicone bordered foam dressing with Super Absorbers			7 cm x 10 cm / oval	7328414	10
<b>Cutimed® Siltec® Sacrum</b>  Self-adherent silicone bordered foam dressing for the sacral area with Super Absorbers			17.5 cm x 17.5 cm	7328700	5
			23 cm x 23 cm	7328701	5
<b>Cutimed® Siltec® Heel</b>  Soft tack silicone foam dressing for the heel with Super Absorbers			16 cm x 24 cm / Heel 2D	7328805	5
			Heel 3D	7328601	5

## BSN medical Australia & New Zealand

L6, 211 Wellington Road,  
Mulgrave, Victoria, Australia 3170

[www.bsnmedical.com.au](http://www.bsnmedical.com.au)  
[www.bsnmedical.co.nz](http://www.bsnmedical.co.nz)

BSN0239

04/2017

Customer Service Australia  
T 1300 276 633 F 1300 998 830

Customer Service New Zealand  
T 0508 276 111 F 0508 998 830

®Registered Trademark

**BSN** medical

an SCA company