

ProShield®

Mask Range

Fluid
Resistant

PROSHIELD® SOFT FR LOOP MASK

FLUID RESISTANT, LEVEL 2,
SUBMICRON FILTRATION,
LOOP MASK

- Fluid resistant
- Ultra soft inner liner suitable for sensitive skin
- Submicron filtration
- Soft ear loops
- Excellent breathability
- Complies with AS4381: 2015

THERAPIES. HAND IN HAND.

BSN medical

an Essity company

PROSHIELD® SOFT FR LOOP MASK

- Submicron filtration
- Fluid Resistance (120mmHg)
- Triple pleated design for easy breathing
- Comfortable hypoallergenic inner liner
- Soft and durable loops

LEVEL 2* PROTECTION

Level 2 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the medium velocity of 120mmHg, bacterial filtration efficiency and differential pressure.

BFE

With a **Bacterial Filtration Efficiency (BFE)** of 98% for media greater than 3 microns this mask provides superior protection from potential infection.

Δ P

The mask design provides excellent breathability and a **low delta P** (≤ 1.5) ensuring coolness and comfort even when worn for an extended period.

PFE

The super high **Particulate Filtration Efficiency (PFE)** material filters >99% of particles greater than 0.1 microns. When securely fitted, the mask helps to provide an effective barrier to particulate matter.



The **Fluid Resistant** qualities of the mask provide protection against fluid (120mmHg/16.0kPa) strikethrough.

APPLICATION

For use in emergency departments, dentistry, changing dressings on small wounds or healing wounds and where minimal blood droplet exposure may possible occur (e.g. endoscopy procedures).



FLUID RESISTANCE

- The fluid resistance of a surgical face mask is the ability of the mask to limit the passage of blood or body fluids from the outer facing through the mask layer to the inner facing.
- The pressure of the spurt is dependent on a number of factors. These include:
 - ▶ Blood pressure of the patient
 - ▶ Velocity of the spurt
 - ▶ Distance between the source of the spurt to the mask wearer

ORDERING INFORMATION

Ref-no.	Description	Items per unit	Shipper Qty
71920-07	ProShield® Soft FR Loop	bx/50	30 bxs

*AS4381:2015 Single-use face masks for use in healthcare

BSN medical Australia & New Zealand

L6, 211 Wellington Road,
Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au
www.bsnmedical.co.nz

BSN0349

10/2017

Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand
T 0508 276 111 F 0508 998 830

BSN medical

an Essity company