

Introducing Tubular - Form SSB™: The ideal shaped tubular compression bandage.







Tubular - FormSSB[™]

Tubular Shaped Support Bandage

Tubular - FormSSB[™] is ideally designed for the treatment of orthopaedic and vascular conditions.

Tubular - FormSSB[™] is used for:

- Acute Venous Insufficiency
- Control of Oedema / Swelling of the limbs
- Joint Effusion and Soft Tissue Injuries
- Varicose Ulcers

Features:

- 100% Australian Made
- Provides Moderate Compression (18-22mmHg) at the ankle and wrist under a single layer
- Product treated with Unique Low Fray Formula
- Easy to use Measuring Guides available
- Unique colour coding system on Full Arm 1/2 Leg sizes

Measuring Guide for Application				
Code	Size	Application	Box	Measurement**
SM700	Small	Full Arm – 1/2 Leg	10	15 - 19cm*
5		3		
SM701	Medium	Full Arm – 1/2 Leg	10	18 - 23cm*
SM702	Large	Full Arm – 1/2 Leg	10	22 - 27cm*
SM707	Extra Large	Full Arm – 1/2 Leg	10	26 - 35cm*
SM708	Extra Extra Large	Full Arm – 1/2 Leg	10	34 - 39cm*
SM703	Small	Full Leg	10	15 - 19cm*
SM704	Medium	Full Leg	10	18 - 23cm*
SM705	Large	Full Leg	10	22 - 27cm*
SM706	Extra Large	Full Leg	10	26 - 35cm*
				LATEX
*When applying to limbs, measure to fit narrowest part.				

Fibre content: Cotton Elastodiene, Polyurethane and Polyester

Sutherland Medical Pty Limited

Unit 2, 195 Chesterville Road Moorabbin Vic 3189

Ph: 1300 664 027

Email: orders@sutherlandmedical.com.au Web Site: www.sutherlandmedical.com.au